

Intelligent Pump UI-22 Series

for continuous flow without pulsation
for Metal-Free application

Ultimate Low Pulsation

Innovative Suction / Discharge Mechanism
without flow restrictions by trapped air bubbles

Advanced Plunger Seal Performance
with introducing high-viscosity liquid such as diesel



Stainless steel spec.



PEEK spec.



PCTFE (Teflon®) spec.

UI-22 series Specifications

Part number	UI-22-110			UI-22-410	
Pump head material	SS	PEEK	PCTFE	SS	PCTFE
Max. pressure	35 MPa (5000 psi)	20 MPa (2900 psi)	2 MPa (290 psi)	5 MPa (725 psi)	2 MPa (290 psi)
Flow rate range	0.001 — 9.999 mL/min			0.01 — 99.9 mL/min	
Communication	RS232 , RS485				

FLOM Corporation

5-32-10 Shinmachi, Ome, Tokyo 198-0024 Japan

Phone: +81-428-30-7454

Website: www.flom.co.jp

Email: info@flom.co.jp